

WEDGES

Wedges are preshaped cancellous or corticocancellous bone grafts, mainly used in corrective osteotomy. C+TBA provides a wide range of wedges (cancellous or cortico-cancellous) with different angles and sizes to precisely address the indication and the patient's individual anatomic preconditions.

Cancellous wedge

Origin: Human
 Tissue: Cancellous bone
 Processing: Allotec® purification procedure
 Inactivation: Min. SAL10⁻⁶ for viruses and bacteria
 Sterilisation: Gamma irradiation
 Application: Corrective osteotomy
 Rehydration: Min. 10 minutes

Cortico-cancellous wedge

Origin: Human
 Tissue: Cortical and cancellous bone
 Processing: Allotec® purification procedure
 Inactivation: Min. SAL10⁻⁶ for viruses and bacteria
 Sterilisation: Gamma irradiation
 Application: Corrective osteotomy
 Rehydration: Min. 10 minutes



Cancellous wedge



ORDER INFORMATION*

DESCRIPTION	A	ITEM NUMBER	S	D [mm]	H [mm]
Cancellous wedge	7°	ALO462	S	<45	5,0
		ALO460	L	≥45	7,0
	10°	ALO465	S	<45	7,0
		ALO463	L	≥45	10,0
	13°	ALO468	S	<45	10,0
		ALO466	L	≥45	13,0
	16°	ALO470	S	<45	13,0
		ALO469	L	≥45	16,0
Cortico-cancellous wedge	15°	ALO410	-	n.a.	10,0



Cortico-cancellous wedge

* Please note: Due to the nature of human bone tissue and the technical possibilities of shaping, slight deviations of the specified sizes may occur.

BT Medicals AB
 nordics & baltics

Scheelevägen 30-32
 223 63 Lund

+46 364404444
 office@btmedical.com

www.btmedical.com